

STOP! READ FIRST!

IMPORTANT: DONORS SHOULD NOT PLACE ANYTHING (INCLUDING FOOD, DRINK, GUM OR TOBACCO PRODUCTS) IN THEIR MOUTH FOR AT LEAST 10 MINUTES PRIOR TO THE PROCEDURE. DO NOT BITE, SUCK, OR CHEW ON THE SPONGE. REFRAIN FROM TALKING WHILE COLLECTION SWAB IS IN THE MOUTH.



Follow these instructions for the best results with this innovative screening device, which includes market-leading drug cutoffs and the first-ever sponge saturation indicator.

1 The purpose of Step 1 is to saturate the sponge with saliva. With collection swab (B), have donor sweep the inside of mouth (cheek, gums, tongue) several times, then hold swab in closed mouth until color on the saturation indicator strip appears in the indicator window. Donor must leave swab in mouth until instructed to remove it. Donor may remove swab from mouth when color appears in the indicator window.

Note: If at 7 minutes, color on saturation indicator has not appeared in indicator window, proceed with test (Step 2).



2 Remove collection swab from mouth and insert sponge first into screening device (A), pushing until the locking flange locks in place in the bottom of device. Once locked in place, the device is airtight, tamper evident, and ready to dispose of after use or send to lab for confirmation (on presumptive positive result). Ensure that specimen is contacting all test strips. If not, rotate the device side to side / front to back to disperse the specimen within the chamber.




*EXAMPLE:
45° maximum
rotation*



With specimen dispersed, set device upright on flat surface. Keep upright while test is running.

3 Read negative results as soon as two lines appear on any test strip, often within 2 minutes. Any indication of a line in the test region (T) should be considered a line, and therefore a negative result. Read presumptive positive results at 10 minutes.




*EXAMPLE:
Presumptive
positive for
cocaine (COC)*



Note: The THC (marijuana) strip has improved sensitivity, which may show a less intense color band in the test region (T). This is an expected and normal response.

Screening Device (A)

Collection Swab (B)



LIMITATIONS OF THE TEST AS PROVIDED BY MANUFACTURER

1. This drug test ("device") is for forensic use, and should be only used for the qualitative detection of drug of abuse in oral fluid.
2. This assay provides a preliminary analytical test result only. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) has been established as the preferred confirmatory method by the National Institute on Drug Abuse (NIDA). Clinical consideration and professional judgment should be applied to any test result, particularly when preliminary positive results are indicated.
3. There is a possibility that technical or procedural errors as well as other substances and factors may interfere with the test and cause false results.
4. A positive result indicates the presence of a drug/metabolite only, and does not indicate or measure intoxication.
5. A negative result does not at any time rule out the presence of drugs/metabolites in urine, as they may be present below the minimum detection level of the test.
6. This test does not distinguish between drugs of abuse and certain medications.

ADDITIONAL LEGAL DISCLAIMER

In making this drug test ("device") available at no cost to you to use at your own direction, Operation UNITE and Appalachia HIDTA, their affiliates, partners, owners, employees, and suppliers (collectively "UNITE and HIDTA") assume no responsibility for consequences resulting from its use. UNITE and HIDTA specifically disclaim any liability, whether based in contract, tort, strict liability or otherwise, for any direct, indirect, incidental, consequential or special damages arising out of in any way connected with access to or use of the device. We assume no responsibility for problems or accidents arising from the use of the device in violation of state laws. No rights or claims can be derived from (the use of) this device. Use of device is at your own risk. A 100% guarantee of accuracy cannot be given nor assured. Actual device packaging and materials may contain more and/or different information than provided by UNITE and HIDTA. We recommend that you do not solely rely on the information presented and that you always read labels, warnings, and directions before using a device. For additional information about the device, please contact the manufacturer. Any content provided by UNITE and HIDTA is for reference purposes and is not intended to substitute for advice given by a physician, pharmacist, or other licensed health-care professional. You should not use this information as self-diagnosis or for treating a health problem or disease. Contact your health-care provider immediately if you suspect that you have a medical problem.